

Rabbit anti-Human CD7(scFv)

Catalog No.: A28936

Basic Information

Observed MW

Calculated MW

25 kDa

Category

Primary antibody

Applications

FC

Cross-Reactivity

Human

CloneNo number

ARC54668 (scFv)

Background

This gene encodes a transmembrane protein which is a member of the immunoglobulin superfamily. This protein is found on thymocytes and mature T cells. It plays an essential role in T-cell interactions and also in T-cell/B-cell interaction during early lymphoid development.

Recommended Dilutions

FC 1:100 - 1:500

Immunogen Information

Gene ID

924

Swiss Prot

P09564

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

GP40; TP41; Tp40; LEU-9

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG


Purification

Affinity purification

Storage

Store at 2-8°C

Validation Data



Flow cytometry: 1×10^6 Human PBMC were surface-stained with Rabbit IgG isotype control (AC042,2 $\mu\text{g/mL}$,left) or Rabbit anti-Human CD7(scFv) mAb (A28936,2 $\mu\text{g/mL}$,right), followed by Alexa Fluor® 647 His tag staining. Cells in the lymphocyte gate were used for analysis.