

Rabbit anti-Human PLXDC1 mAb

Catalog No.: A28975

Basic Information

Observed MW

Calculated MW

Category

Primary antibody

Applications

FC

Cross-Reactivity

Human

CloneNo number

ARC5246

Background

Recommended Dilutions

FC 1:100-1:500

Immunogen Information

Gene ID

57125

Swiss Prot

Q8IUK5

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: 25mM HEPES, 150mM NaCl, 0.03% Proclin, 50% glycerol, pH8.0.