

APC Mouse IgG 1 isotype control

Catalog No.: A29019

Basic Information

Observed MW

Calculated MW

Category

Primary antibody

Applications

FC

Cross-Reactivity

CloneNo number

AMC50005

Conjugate

APC. Ex:650nm. Em:660nm.

Background

Recommended Dilutions

FC ≤0.25 µg per million cells in
100 µl volume

Immunogen Information

Gene ID

Swiss Prot

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Mouse

Isotype

IgG


Purification

Affinity purification

Storage

Store at -20°C

Validation Data



Flow cytometry: 1X10⁶ Human PBMC were surface-stained with APC Mouse IgG1 isotype control (A29019,0.25 µg,left) or APC Mouse anti-Human TCR α/β mAb (A27980,5 µl/Test,right). Cells in the lymphocyte gate were used for analysis.