

# PE-ABflo® 594 anti-Mouse CD11c mAb

Catalog No.: A29386

## Basic Information

### Observed MW

### Calculated MW

### Category

Primary antibody

### Applications

FC

### Cross-Reactivity

Mouse

### CloneNo number

ACC10002

### Conjugate

PE-ABflo® 594. Ex:566nm. Em:615nm.

## Background

Predicted to enable receptor tyrosine kinase binding activity. Involved in positive regulation of gene expression and positive regulation of myelination. Acts upstream of or within defense response to virus. Located in external side of plasma membrane. Is expressed in heart and hemolymphoid system. Orthologous to human ITGAX (integrin subunit alpha X).

## Recommended Dilutions

**FC** ≤0.25 µg per million cells in  
100 µl volume

## Immunogen Information

### Gene ID

16411

### Swiss Prot

Q9QXH4

### Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

### Synonyms

Cr4; N418; Cd11c

## Contact

☎ | 400-999-6126

✉ | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

🌐 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Product Information

### Source

Armenian Hamster

### Isotype

Armenian hamster / IgG

### Purification


Affinity purification

### Storage

Store at 2-8°C. Avoid freeze.

Buffer: PBS with 0.09% Sodium azide, 0.2% BSA, pH7.3.

# Validation Data



Flow cytometry: 1X10<sup>6</sup> C57BL/6 mouse splenocytes were surface-stained with ABflo® 647 Rat anti-Mouse MHC Class II (I-A/I-E) mAb (A25482,5 µl/Test) and PE-ABflo® 594 Armenian Hamster IgG isotype control (5 µl/Test,left) or PE-ABflo® 594 anti-Mouse CD11c mAb (A29386,0.25 µg,right).