

ABflo® 647 Rabbit anti-Mouse CD366/TIM-3/HAVCR2 mAb

Catalog No.: A29521

Basic Information

Observed MW

Calculated MW

Category

Primary antibody

Applications

FC

Cross-Reactivity

Mouse

CloneNo number

ARC59425

Conjugate

ABflo® 647. Ex:648nm. Em:664nm.

Background

Predicted to enable metal ion binding activity. Involved in several processes, including regulation of cytokine production; regulation of leukocyte activation; and toll-like receptor signaling pathway. Located in cell surface; early endosome; and immunological synapse. Is expressed in several structures, including dorsal aorta; gut; liver; lung; and reproductive system. Human ortholog(s) of this gene implicated in hepatocellular carcinoma and renal cell carcinoma. Orthologous to human HAVCR2 (hepatitis A virus cellular receptor 2).

Recommended Dilutions

FC ≤0.5 µg per million cells in
100 µl volume

Immunogen Information

Gene ID

171285

Swiss Prot

Q8VIM0

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

Tim3; TIM-3; Timd3

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at 2-8°C. Avoid freeze.

Buffer: PBS with 0.09% Sodium azide, 0.2% BSA, pH7.3.