

RelB Rabbit mAb

Catalog No.: A29551 **Recombinant**

Basic Information

Observed MW

70kDa

Calculated MW

62 kDa

Category

Primary antibody

Applications

WB,IP,IF-F,IF-P,IHC-P,ELISA

Cross-Reactivity

Human, Mouse, Rat

CloneNo number

ARC4182

Background

Enables RNA polymerase II cis-regulatory region sequence-specific DNA binding activity and protein kinase binding activity. Involved in lymphocyte differentiation and negative regulation of interferon-beta production. Located in cytosol and nucleoplasm. Part of chromatin; nucleus; and transcription repressor complex. Colocalizes with centrosome. Implicated in breast cancer and immunodeficiency 53. Biomarker of breast cancer and transitional cell carcinoma.

Recommended Dilutions

WB 1:50000-1:200000

IP 0.5µg-4µg antibody for
200µg-400µg extracts of
whole cells

IF-F 1:1000-1:2000

IF-P 1:1000-1:2000

IHC-P 1:8000-1:32000

ELISA Recommended starting
concentration is 1 µg/mL.
Please optimize the
concentration based on your
specific assay requirements.

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

Immunogen Information

Gene ID

5971

Swiss Prot

Q01201

Immunogen

This information is considered to be commercially sensitive.

Synonyms

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C



www.abclonal.com.cn
