

DYNC1I1 Rabbit pAb

Catalog No.: A3833

Basic Information

Observed MW

73kDa

Calculated MW

73kDa

Category

Primary antibody

Applications

WB, IHC-P, ELISA

Cross-Reactivity

Human, Mouse, Rat

Background

Enables spectrin binding activity. Involved in vesicle transport along microtubule. Located in several cellular components, including kinetochore; recycling endosome; and spindle pole.

Recommended Dilutions

WB 1:200 - 1:2000

IHC-P 1:20 - 1:200

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

1780

Swiss Prot

O14576

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

DNCI1; DNCIC1; DYNC1I1

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG



Purification

Affinity purification


Storage

Store at -20°C

Validation Data



Immunohistochemistry analysis of paraffin-embedded Human stomach using DYNC111 Rabbit pAb (A3833) at dilution of 1:100 (40x lens). Microwave antigen retrieval performed with 0.01M PBS Buffer (pH 7.2) prior to IHC staining.



Immunohistochemistry analysis of paraffin-embedded Rat brain using DYNC111 Rabbit pAb (A3833) at dilution of 1:100 (40x lens). Microwave antigen retrieval performed with 0.01M PBS Buffer (pH 7.2) prior to IHC staining.