

CNNM3 Rabbit pAb

Catalog No.: A4625

Basic Information

Observed MW

Refer to figures

Calculated MW

76kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

Predicted to enable transmembrane transporter activity. Predicted to be involved in ion transport; magnesium ion homeostasis; and transmembrane transport. Located in membrane.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

26505

Swiss Prot

Q8NE01

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

ACDP3; CNNM3

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.