

Rat Control IgG

Catalog No.: AC034

Basic Information

Observed MW

Calculated MW

Category

Loading control antibody

Applications

IP

Cross-Reactivity

-

Background

Recommended Dilutions

IP 1:50 - 1:200

Immunogen Information

Gene ID

Swiss Prot

Immunogen

This information is considered to be commercially sensitive.

Synonyms

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rat

Isotype

IgG

Purification

Potein A/G purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.