

ABflo® 594 Rabbit anti His-Tag mAb

Catalog No.: AE152

Basic Information

Observed MW

Calculated MW

Category

Tag antibody

Applications

FC (intra),FC

Cross-Reactivity

Species independent

CloneNo number

ARC2696RE

Conjugate

ABflo® 594. Ex:588nm. Em:604nm.

Background

Consecutive histidine residues (usually 6 to 10 in length) are often inserted into the amino acid sequences of recombinant proteins. The resulting His-tagged proteins can be detected by using His antibodies.

Recommended Dilutions

FC (intra) FC≤0.25 µg per million
cells in 100 µl volume

Immunogen Information

Gene ID

Swiss Prot

Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

Synonyms

His;His tag;His-tag

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at 2-8°C. Avoid freeze.

Buffer: PBS with 0.09% Sodium azide,0.2% BSA,pH7.3.