

Phospho-IKK β -Y199 Rabbit pAb

Catalog No.: AP0369

1 Publications

Basic Information

Observed MW

Calculated MW

87kDa

Category

Primary antibody

Applications

Cross-Reactivity

Human, Mouse, Rat

Background

The protein encoded by this gene phosphorylates the inhibitor in the inhibitor/NF-kappa-B complex, causing dissociation of the inhibitor and activation of NF-kappa-B. The encoded protein itself is found in a complex of proteins. Several transcript variants, some protein-coding and some not, have been found for this gene.

Recommended Dilutions

Immunogen Information

Gene ID

3551

Swiss Prot

O14920

Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

Synonyms

IKK2; IKKB; IMD15; IMD15A; IMD15B; NFKB1KB; IKK-beta; Phospho-IKK β -Y199

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.