

# ABflo® 647-conjugated Rabbit anti-Rat IgG (H+L)

Catalog No.: AS068 2 Publications

## **Basic Information**

#### **Observed MW**

**Calculated MW** 

## Category

Secondary antibody

## **Applications**

IF/ICC,FC

## **Cross-Reactivity**

Rat

#### Conjugate

ABflo® 647. Ex:648nm. Em:664nm.

# **Background**

Secondary antibodies are affinity-purified antibodies which will work with target-specific primary antibody in the detection, sorting or purification of its specified target. Secondary antibodies offer increased versatility enabling users to use many detection systems (e.g. HRP, AP, fluorescence). They can also provide greater sensitivity through signal amplification as multiple secondary antibodies . Most commonly, secondary antibodies are generated by immunizing the host animal (different from host species of primary antibody) with a pooled population of normal immunoglobulins from the host species of primary antibody and can be further purified and modified (i.e. antibody fragmentation, label conjugation, etc.) to ensure well-characterized specificity to corresponding normal immunoglobulins.

## **Recommended Dilutions**

**IF/ICC** 1:200 - 1:500

FC 1:500 - 1:2000

# **Immunogen Information**

Gene ID Swiss Prot

#### **Immunogen**

This information is considered to be commercially sensitive.

#### **Synonyms**

## **Contact**

<b>a</b>	400-999-6126
×	cn.market@abclonal.com.cn
	www.abclonal.com.cn

# **Product Information**

SourceIsotypePurificationRabbitIgGAffinity purification

#### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS containing 50% glycerol, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.

# **Validation Data**





Flow cytometry: 1X10^6 RK13 cells (negative control,left) and RK13-CD20 transfection cells (right) were surface-stained with rat anti-mouse CD20 Antibody (1:100,orange line) or secondary antibody only (blue line). Non-fluorescently stained RK13 and RK13 transfection cells were used as blank control (red line). ABflo® 647-conjugated Rabbit Anti-Rat IgG (H+L)(AS068, 1:1000) was used as a secondary antibody.