

HRP-conjugated Rabbit anti-Human IgM Rabbit mAb

Catalog No.: AS110 **Recombinant**

Basic Information

Observed MW

Refer to figures

Calculated MW

Category

Secondary antibody

Applications

ELISA

Cross-Reactivity

CloneNo number

ARC63180-02-HRP

Conjugate

HRP

Background

Recommended Dilutions

ELISA 1:5000-1:10000

Immunogen Information

Gene ID

Swiss Prot

Immunogen

Synonyms

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype


Purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.

Affinity of Rabbit anti Rabbit Anti-Human IgM mAb, HRP conjugated



Dose response curve of HRP conjugated Human IgM mAb measured by ELISA. 1 µg/mL of various immunoglobulins were coated to 384-well plate., blank wells without protein were used as negative control (NC). The coated plate was blocked and subsequently incubated with 25 µL of HRP conjugated Rabbit Anti-Human IgM mAb in a 3 fold serial dilution from 1 µg/mL to 1.52