

## Streptavidin-HRP Diluent (R3)

**Catalog :** RM00025

**Product Description:** Diluent for Streptavidin-HRP

**Storage:**Store at 2 - 8°C

### **Handling:**

- 1.Bring all reagents to room temperature before use.
- 2.Refer to the Product Information Sheet for the applicable ELISA kit.

### **PRECAUTIONS:**

- 1.Laboratory supplies should be specific to avoid cross contamination.
- 2.Wear protective gloves, clothing, eye, and face protection. Wash hands thoroughly after handling.

**\*FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.**