

Rabbit anti-Human Galectin 1/LGALS1 mAb (CAP)

Catalog No.: RM17806

Basic Information

Catalog No.
RM17806

Catagory
Elisa Antibody Kit

Application
ELISA

Product Information

Ig Type
Rabbit IgG

Purification
Affinity purification

Endotoxin Level

Storage
Store at -20°C. Avoid freeze / thaw cycles. Preservative 0.05%ProClin 300.
Avoid repeated freeze-thaw cycles.

Formulation
Supplied as a 0.2um filtered solution in PBS with 0.05%ProClin 300,PH 7.4.

Contact

 | order@abclonal.com

 | support@abclonal.com

 | www.abclonal.com

Background

The galectins are a family of beta-galactoside-binding proteins implicated in modulating cell-cell and cell-matrix interactions. This gene product may act as an autocrine negative growth factor that regulates cell proliferation.

Immunogen Information

Immunogen
Recombinant Human Galectin 1/LGALS1

Cross-Reactivity


/

Assay Applications


Human Galectin 1/LGALS1 Sandwich ELISA Immunoassay

	Recommended Concentration	Sample
ELISA Capture	0.25-1ug/mL	Rabbit anti-Human Galectin 1/LGALS1 (CAP)(Cat. No.RM17806)
ELISA Detection	0.0625-0.25ug/mL	Rabbit anti-Human Galectin 1/LGALS1 (DET)(Cat. No.RM17807)
Standard	0.78-50ng/mL	Recombinant Human Galectin 1/LGALS1 Protein(Cat. No.RP00007)

Validation Data



This standard curve is only for demonstration purposes. A standard curve should be generated for each assay.



Placed at 37°C for 3 days, the stability of the standard curve all conform to CV <10%.