

Rabbit anti-Human PGI₂Ab (CAP)

Catalog No.: RM17939

Basic Information

Catalog No.
RM17939

Catagory
Elisa Antibody Kit

Application
ELISA

Product Information

Ig Type
Rabbit IgG

Purification
Affinity purification

Endotoxin Level

Storage
Store at -20°C. Avoid freeze / thaw cycles. Preservative 0.05%ProClin 300.
Avoid repeated freeze-thaw cycles.

Formulation
Supplied as a 0.2um filtered solution in PBS with 0.05%ProClin 300,PH 7.4.

Contact

 | order@abclonal.com

 | support@abclonal.com

 | www.abclonal.com

Background

This gene encodes a protein precursor of the digestive enzyme pepsin, a member of the peptidase A1 family of endopeptidases. The encoded precursor is secreted by gastric chief cells and undergoes autocatalytic cleavage in acidic conditions to form the active enzyme, which functions in the digestion of dietary proteins. This gene is found in a cluster of related genes on chromosome 11, each of which encodes one of multiple pepsinogens. Pepsinogen levels in serum may serve as a biomarker for atrophic gastritis and gastric cancer.

Immunogen Information

Immunogen
Recombinant Human PGI

Cross-Reactivity


/

Assay Applications

Human PGI Sandwich ELISA Immunoassay

	Recommended Concentration	Sample
ELISA Capture	1-4ug/mL	Rabbit anti-Human PGI (CAP)(Cat. No.RM17939)
ELISA Detection	0.17-0.042ug/mL	Rabbit anti-Human PGI (DET)(Cat. No.RM17940)
Standard	3.91-250pg/mL	Recombinant Human PGI Protein

Validation Data



This standard curve is only for demonstration purposes. A standard curve should be generated for each assay.