

Rabbit anti-Human SIGLEC5 mAb(CAP)

Catalog No.: RM18041

Basic Information

Catalog No.
RM18041

Catagory
Elisa Antibody Kit

Application
ELISA

Product Information

Ig Type
Rabbit IgG

Purification
Affinity purification

Endotoxin Level

Storage
Store at -20°C. Avoid freeze / thaw cycles. Preservative 0.05%ProClin 300.
Avoid repeated freeze-thaw cycles.

Formulation
Supplied as a 0.2umfiltered solution in PBS with0.05%ProClin 300 ,PH 7.4.

Contact

 | order@abclonal.com

 | support@abclonal.com

 | www.abclonal.com

Background

This gene encodes a member of the sialic acid-binding immunoglobulin-like lectin (Siglec) family. These cell surface lectins are characterized by structural motifs in the immunoglobulin (Ig)-like domains and sialic acid recognition sites in the first Ig V set domain. The encoded protein is a member of the CD33-related subset of Siglecs and inhibits the activation of several cell types including monocytes, macrophages and neutrophils. Binding of group B Streptococcus (GBS) to the encoded protein plays a role in GBS immune evasion.

Immunogen Information

Immunogen
Recombinant Human SIGLEC5

Cross-Reactivity


/

Assay Applications

Human SIGLEC5 Sandwich ELISA Immunoassay

	Recommended Concentration	Sample
ELISA Capture	1-4ug/mL	Rabbit anti-Human SIGLEC5 mAb(CAP)(Cat. No.RM18041)
ELISA Detection	0.025-0.1ug/mL	Rabbit anti-Human SIGLEC5 mAb(DET)(Cat. No.RM18042)
Standard	46.88-3000pg/mL	Recombinant Human SIGLEC5 Protein

Validation Data



This standard curve is only for demonstration purposes. A standard curve should be generated for each assay.