

Sheep anti-General Estradiol mAb

Catalog No.: RMK0469

Basic Information

Catalog No.

RMK0469

Category

Elisa Antibody Kit

Application

ELISA

Product Information

Ig Type

Sheep IgG

Purification

Affinity purification

Endotoxin Level

Storage

Store at -20°C

Avoid repeated freeze-thaw cycles.

Formulation

Supplied as a 0.2um filtered solution in PBS with 0.05%ProClin 300,PH 7.4.

Contact

 | order@abclonal.com

 | support@abclonal.com

 | www.abclonal.com

Background

Immunogen Information

Immunogen

Recombinant Estradiol protein

Cross-Reactivity


/

Assay Applications

General Estradiol Competitive ELISA

		Recommended Concentration	Sample
ELISA	Capture	0.25-1 ug/mL	Recombinant General Estradiol Protein
	Detection	1-5 ng/mL	Sheep anti-General Estradiol mAb (Cat. No. RMK0469)
	Standard	0.014-10 ng/mL	Recombinant General Estradiol Protein

Validation Data



This standard curve is only for demonstration purposes. A standard curve should be generated for each assay.