# Mouse anti-Human IgA mAb



Catalog No.: RMK0574

## **Basic Information**

## Catalog No.

RMK0574

#### Catagory

Elisa Antibody Kit

# **Application**

**ELISA** 

## **Product Information**

## **Ig Type**

Mouse IgG2b

## **Purification**

Affinity purification

### **Endotoxin Level**

#### Storage

Store at -20°C.

Avoid repeated freeze-thaw cycles.

## **Formulation**

Supplied as a 0.2  $\mu m$  filtered solution in PBS with 0.05% Proclin 300, PH 7.4.

## Contact

•	order@abclonal.com
<u>a</u>	support@abclonal.com
•	www.abclonal.com

# **Background**

Contributes to immunoglobulin receptor binding activity. Involved in antibacterial humoral response; glomerular filtration; and positive regulation of respiratory burst. Located in extracellular space. Part of monomeric IgA immunoglobulin complex and secretory dimeric IgA immunoglobulin complex.

# **Immunogen Information**

## **Immunogen**

## **Cross-Reactivity**

\_ \_ \_

# **Assay Applications**