

Mouse anti-Human IgA mAb

Catalog No.: RMK0575

Basic Information

Catalog No.

RMK0575

Catagory

Elisa Antibody Kit

Application

ELISA

Product Information

Ig Type

Mouse IgG1

Purification

Affinity purification

Endotoxin Level

Storage

Store at -20°C.

Avoid repeated freeze-thaw cycles.

Formulation

Supplied as a 0.2 μm filtered solution in PBS with 0.05% Proclin 300, PH 7.4.

Contact

<u>.</u>	order@abclonal.com
<u>a</u>	support@abclonal.com
•	www.abclonal.com

Background

Contributes to immunoglobulin receptor binding activity. Involved in antibacterial humoral response; glomerular filtration; and positive regulation of respiratory burst. Located in extracellular space. Part of monomeric IgA immunoglobulin complex and secretory dimeric IgA immunoglobulin complex.

Immunogen Information

Immunogen

Cross-Reactivity

Assay Applications