

Mouse anti-Human IgE mAb

Catalog No.: RMK0576

Basic Information

Catalog No.
RMK0576

Category
Elisa Antibody Kit

Application
ELISA

Product Information

Ig Type
Mouse IgG2a


Purification
Affinity purification


Endotoxin Level

Storage
Store at -20°C.
Avoid repeated freeze-thaw cycles.

Formulation
PBS, pH 7.4, 0.09 % sodium azide

Contact

 | order@abclonal.com

 | support@abclonal.com

 | www.abclonal.com

Background

IgE is the class of antibodies produced in the lungs, skin, and mucous membranes. It may protect against parasite invasion, but it is a major factor in allergic reactions. The antigen-specific IgE interacts with mast cells and eosinophils, triggers the release of histamine, leukotrienes and other substances that lead to the itching, sneezing and congestion of allergies - and the life threatening respiratory distress of asthma and anaphylactic shock.

Immunogen Information

Immunogen

Cross-Reactivity

/

Assay Applications
