

Rabbit anti-Human CCL20/MIP-3 alpha mAb

Catalog No.: RMK0592

Basic Information

Catalog No.
RMK0592

Catagory
Elisa Antibody Kit

Application

Product Information

Ig Type
Rabbit IgG


Purification
Affinity purification


Endotoxin Level
/


Storage
Store at -20°C.
Avoid repeated freeze-thaw cycles.

Formulation
Supplied as a 0.2µm filtered solution in PBS with 0.05% Proclin300, PH 7.4.

Contact

 | order@abclonal.com

 | support@abclonal.com

 | www.abclonal.com

Background

This antimicrobial gene belongs to the subfamily of small cytokine CC genes. Cytokines are a family of secreted proteins involved in immunoregulatory and inflammatory processes. The CC cytokines are proteins characterized by two adjacent cysteines. The protein encoded by this gene displays chemotactic activity for lymphocytes and can repress proliferation of myeloid progenitors. Two transcript variants encoding different isoforms have been found for this gene.

Immunogen Information

Immunogen
Recombinant Human CCL20/MIP-3 alpha protein

Cross-Reactivity

Assay Applications