

Rabbit anti-Human CCL20/MIP-3 alpha mAb

Catalog No.: RMK0596

Basic Information

Catalog No.
RMK0596

Catagory
Elisa Antibody Kit

Application

Product Information

Ig Type
Rabbit IgG

Purification
Affinity purification

Endotoxin Level
/

Storage
Store at -20°C.
Avoid repeated freeze-thaw cycles.

Formulation
Supplied as a 0.2µm filtered solution in PBS with 0.05% Proclin300, PH 7.4.

Contact

		order@abclonal.com
		support@abclonal.com
		www.abclonal.com

Background

This antimicrobial gene belongs to the subfamily of small cytokine CC genes. Cytokines are a family of secreted proteins involved in immunoregulatory and inflammatory processes. The CC cytokines are proteins characterized by two adjacent cysteines. The protein encoded by this gene displays chemotactic activity for lymphocytes and can repress proliferation of myeloid progenitors. Two transcript variants encoding different isoforms have been found for this gene.

Immunogen Information

Immunogen
Recombinant Human CCL20/MIP-3 alpha protein

Cross-Reactivity

Assay Applications

Human CCL20/MIP-3 alpha Sandwich Immunoassay

		Recommended Concentration	Sample
Multiplex	Capture	3-20 µg/mL	Rabbit anti-Human CCL20/MIP-3 alpha mAb (Cat. No. RMK0592)
	Detection	0.017-2 µg/mL	Rabbit anti-Human CCL20/MIP-3 alpha mAb (Cat. No. RMK0596)
	Standard	6.86-5000 pg/mL	Recombinant Human CCL20/MIP-3 alpha protein