

Rabbit anti-Mouse CXCL10/IP-10 mAb

Catalog No.: RMK0608

Basic Information

Catalog No.
RMK0608

Catagory
Elisa Antibody Kit

Application

Product Information

Ig Type
Rabbit IgG

Purification
Affinity purification

Endotoxin Level
N/A

Storage
Store at -20°C.
Avoid repeated freeze-thaw cycles.

Formulation
Supplied as a 0.2 µm filtered solution in PBS with 0.05% Proclin 300, pH 7.3.

Contact

		order@abclonal.com
		support@abclonal.com
		www.abclonal.com

Background

Predicted to enable chemoattractant activity; chemokine receptor binding activity; and heparin binding activity. Involved in cellular response to interleukin-17. Acts upstream of or within several processes, including defense response to virus; negative regulation of myoblast differentiation; and negative regulation of myoblast fusion. Located in external side of plasma membrane and extracellular space. Is expressed in several structures, including alimentary system; axial skeleton; hemolymphoid system; lacrimal gland; and liver. Human ortholog(s) of this gene implicated in hepatitis B; middle cerebral artery infarction; and type 1 diabetes mellitus. Orthologous to human CXCL10 (C-X-C motif chemokine ligand 10).

Immunogen Information

Immunogen
Recombinant Mouse CXCL10 protein

Cross-Reactivity

Assay Applications

Mouse CXCL10/IP10 Sandwich ELISA

	Recommended Concentration	Sample
ELISA Capture	1-4 µg/mL	Rabbit anti-Mouse CXCL10/IP10 mAb (Cat. No. RMK0608)
ELISA Detection	0.005-0.02 µg/mL	Rabbit anti-Mouse CXCL10/IP10 mAb (Cat. No. RMK0606)
Standard	1.56-100 pg/mL	Recombinant Mouse CXCL10/IP10 Protein

