

Rabbit anti-Mouse Leptin mAb

Catalog No.: RMK0628

Basic Information

Catalog No.
RMK0628

Catagory
Elisa Antibody Kit

Application
Multiplex

Product Information

Ig Type
Rabbit IgG

Purification
Affinity purification

Endotoxin Level
N/A

Storage
Store at -20°C.
Avoid repeated freeze-thaw cycles.

Formulation
Supplied as a 0.2 µm filtered solution in PBS with 0.05% Proclin 300, pH 7.3.

Contact

 | order@abclonal.com

 | support@abclonal.com

 | www.abclonal.com

Background

Immunogen Information

Immunogen
Recombinant Mouse Leptin protein

Cross-Reactivity

Mouse

Assay Applications

Mouse Leptin Sandwich Immunoassay

		Recommended Concentration	Sample
Multiplex	Capture	3-20 µg/mL	Rabbit anti-Mouse Leptin mAb (Cat. No. RMK0625)
	Detection	0.017-2 µg/mL	Rabbit anti-Mouse Leptin mAb (Cat. No. RMK0626)
	Standard	6.86-5000 pg/mL	Recombinant Mouse Leptin protein