

Recombinant Human FKBP1A/FKBP12 Protein

Catalog No.: RP00051 Recombinant

Sequence Information

Species Gene ID Swiss Prot Human 2280 P62942

Tags

No tag

Synonyms

FKBP-12; FKBP-1A; FKBP1; FKBP12; PKC12; PKCI2; PPIASE; FKBP1A; FKBP-1A; FKBP1; FKBP12; P KC12; PKC12; PPIASE

Product Information

Source Purification
E. coli ≥ 90 % as
determined by

determined by SDS-PAGE.

Calculated MW Observed MW

11.95 kDa 14 kDa

Endotoxin

< 0.1 EU/ μ g of the protein by LAL method.

Formulation

Lyophilized from a 0.22 µm filtered solution of 20mM Tris, 150mM NaCl, pH 8.0.Contact us for customized product form or formulation.

Reconstitution

Centrifuge the vial before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stablizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

Contact

Background

The protein is a member of the immunophilin protein family, which play a role in immunoregulation and basic cellular processes involving protein folding and trafficking. The protein is a cis-trans prolyl isomerase that binds the immunosuppressants FK506 and rapamycin. It interacts with several intracellular signal transduction proteins including type I TGF-beta receptor. It also interacts with multiple intracellular calcium release channels, and coordinates multi-protein complex formation of the tetrameric skeletal muscle ryanodine receptor.

Basic Information

Description

Recombinant Human FKBP1A/FKBP12 Protein is produced by *E. coli* expression system. The target protein is expressed with sequence (Met1-Glu108) of human FKBP12 (Accession #NP_000792.1).

Bio-Activity

Shipping

The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below.

Storage

Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.


After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.

Avoid repeated freeze/thaw cycles.

Operational Notes

For your safety and health, please wear a lab coat and disposable gloves for handling.

Validation Data



Recombinant Human FKBP1A/FKBP12 Protein was determined by SDS-PAGE under reducing conditions with Coomassie Blue.