Recombinant Human ALK-4/ACVR1B Protein

ABclonal

www.abclonal.com

Catalog No.: RP00074 Recombinant

Sequence Information

Species Gene ID **Swiss Prot** 91 Human P36896-1

Tags

C-hFc[]his

Synonyms

ACTRIB; ACVRLK4; ALK4; SKR2;ACVR1B;ACVRLK4;ALK4;SKR2

Product Information

Purification HEK293 cells ≥ 95 % as

determined by SDS-PAGE.

Calculated MW Observed MW

40.37 kDa 45-50 kDa

Endotoxin

< 0.1 EU/µg of the protein by LAL method.

Formulation

Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4.Contact us for customized product form or formulation.

Reconstitution

Centrifuge the vial before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stablizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

Contact

<u>a</u>		400-999-6126
\bowtie		cn.market@abclonal.com.cn
•	T	www.abclonal.com.cn

Background

ALK-4 (Activin Receptor-Like Kinase 4) or ACVR1B (Activin A Receptor, type 1B), belongs to the protein kinase superfamily, TKL Ser/Thr protein kinase family, and TGFB receptor subfamily. ALK-4/ACVR1B acts as a transducer of activin or activin like ligands signals. Activin binds to either ACVR2A or ACVR2B and then forms a complex with ACVR1B. The known type II activin receptors include ActRII and ActRIIB, while the main type I activin receptor in mammalian cells is ALK-4 (ActRIB). In the presence of activin, type II and type I receptors form complexes whereby the type II receptors activate ALK-4 through phosphorylation. The activated ALK-4, in turn, transduces signals downstream by phosphorylation of its effectors, such as Smads, to regulate gene expression and affect cellular phenotype. ALK-4/ACVR1B is an important regulator of vertebrate development, with roles in mesoderm induction, primitive streak formation, gastrulation, dorsoanterior patterning, and left-right axis determination.

Basic Information

Description

Recombinant Human ALK-4/ACVR1B Protein is produced by HEK293 cells expression system. The target protein is expressed with sequence (Met 1-Glu 126) of human ACVR1B (Accession #NP_004293.1) fused with a C-hFc his at the C-terminus.

Bio-Activity

Measured by its binding ability in a functional ELISA. Immobilized Human ACVR1B at 0.5 μg/mL (100 μL/well) can bind Human ACVR2B with a linear range of 2.0-286.1

Shipping

The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below.

Storage

Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.


After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.

Avoid repeated freeze/thaw cycles.


Operational Notes

For your safety and health, please wear a lab coat and disposable gloves for handling.

Validation Data



Recombinant Human ALK-4/ACVR1B Protein was determined by SDS-PAGE under reducing conditions with Coomassie Blue.



Immobilized Human ACVR1B at $0.5\mu g/mL$ (100 $\mu L/well$) can bind Human ACVR2B with a linear range of 2.0-286.1ng/mL.