

# Biotinylated Recombinant Human LILRB2/ILT-4/CD85d Protein

Catalog No.: RP00341B **Recombinant**

## Sequence Information

Species	Gene ID	Swiss Prot
Human	10288	Q8N423-1

### Tags

C-His&Avi

### Synonyms

CD85d; ILT4; ILT-4; ILT4CD85d; LILRB2; LIR2; MIR10

## Product Information

### Source

HEK293 cells

### Purification

≥ 95 % as determined by SDS-PAGE; ≥ 95 % as determined by HPLC.

### Calculated MW

50.2 kDa

### Observed MW

70-75 kDa

### Endotoxin

< 1 EU/μg of the protein by LAL method.

### Formulation

Lyophilized from 0.22 μm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.

### Reconstitution

Centrifuge the vial before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

## Contact

☎ | 400-999-6126

✉ | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

## Background

The immunoglobulin-like transcript (ILT) comprise a family of activating and inhibitory type immunoreceptors whose genes are located in the same locus that encodes killer cell Ig-like receptors (KIR). ILT4, also known as LIR-2 and LILRB2, is a type I transmembrane protein expressed primarily on monocytes and dendritic cells (DC). LILRB2 is a receptor for class I MHC antigens. Recognizes a broad spectrum of HLA-A, HLA-B, HLA-C, HLA-G and HLA-F alleles.

## Basic Information

### Description

Biotinylated Recombinant Human LILRB2/ILT-4/CD85d Protein is produced by HEK293 cells expression system. The target protein is expressed with sequence (Gln22-His458) of Human LILRB2/CD85d/ILT4 (Accession #Q8N423-1) fused with a C-His&Avi tag at the C-terminus.

### Bio-Activity

### Shipping

The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below.

### Storage

Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.

After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.


Avoid repeated freeze/thaw cycles.

### Operational Notes


For your safety and health, please wear a lab coat and disposable gloves for handling.




## Validation Data



Biotinylated Recombinant Human LILRB2/ILT-4/CD85d Protein was determined by Tris-Bis PAGE under reducing conditions.



The purity of Biotinylated Human LILRB2 is greater than 95% as determined by SEC-HPLC.



Immobilized Biotinylated Human LILRB2, His Tag at 0.5 µg/mL (100 µL/Well) on streptavidin (5 µg/mL) precoated plate. Dose response curve for Anti-LILRB2 Antibody, hFc Tag with the EC<sub>50</sub> of 4.0ng/mL determined by ELISA.