

# Recombinant Rat TNF-alpha/TNF $\alpha$ /TNFSF2 Protein

Catalog No.: RP00748 **Recombinant**

## Sequence Information

Species	Gene ID	Swiss Prot
Rat		P16599

### Tags

No tag

### Synonyms

Tumor Necrosis Factor; Cachectin; TNF-Alpha; Tumor Necrosis Factor Ligand Superfamily Member2; TNF-a; Tumor Necrosis Factor; Membrane Form; Tumor Necrosis Factor; Soluble Form; Tnf;Tnfa; Tnfsf2

## Product Information

Source	Purification
<i>E. coli</i>	> 95% by SDS-PAGE.

### Endotoxin

< 1 EU/ $\mu$ g of the protein by LAL method.

### Formulation

Lyophilized from a 0.2  $\mu$ m filtered solution of PBS, pH7.4. Contact us for customized product form or formulation.

### Reconstitution

Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water.

## Background

Tumor necrosis factor alpha (TNF-alpha, TNFSF2) is the prototypic ligand of the TNF superfamily. Rat TNF-alpha consists of a 35 amino acid (aa) cytoplasmic domain, a 21 aa transmembrane segment, and a 179 aa extracellular domain (ECD). Within the ECD, rat TNF-alpha shares 94% aa sequence identity with mouse. TNF-alpha is produced by a wide variety of immune, epithelial, endothelial, and tumor cells. TNF exists as a homotrimer and interacts with SPPL2B. TNF is mainly secreted by macrophages and can induce cell death of certain tumor cell lines. TNF is a key cytokine in the development of several inflammatory disorders. It contributes to the development of type 2 diabetes through its effects on insulin resistance and fatty acid metabolism.

## Basic Information

### Description

Recombinant Rat TNF-alpha/TNF $\alpha$ /TNFSF2 Protein is produced by *E. coli* expression system. The target protein is expressed with sequence (Leu80-Leu235) of rat TNF-alpha/TNF $\alpha$ /TNFSF2 (Accession #P16599) fused with an initial Met at the N-terminus.

### Bio-Activity

### Storage

Store the lyophilized protein at -20°C to -80 °C for long term.

After reconstitution, the protein solution is stable at -20 °C for 3 months, at 2-8 °C for up to 1 week.

Avoid repeated freeze/thaw cycles.

## Contact


☎ | 400-999-6126

✉ | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

🌐 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Validation Data

---



Recombinant Rat TNF-alpha/TNF $\alpha$ /TNFSF2 Protein was determined by SDS-PAGE under reducing conditions with Coomassie Blue.