

Recombinant Macaca mulatta SEMG2 Protein

Catalog No.: RP00786

Recombinant

Sequence Information

Species	Gene ID	Swiss Prot
Macaca mulatta		Q95196

Tags
C-His

Synonyms
SEMG2, Semenogelin-2, Semenogelin II, SGII

Product Information

Source	Purification
HEK293 cells	

Calculated MW **Observed MW**

Endotoxin

Formulation

Reconstitution

Background

Participates in the formation of a gel matrix (sperm coagulum) entrapping the accessory gland secretions and ejaculated spermatozoa.

Basic Information

Description

Recombinant Macaca mulatta SEMG2 Protein is produced by HEK293 cells expression system. The target protein is expressed with sequence (Q24-T706) of Macaca mulatta SEMG2 (Accession #Q95196) fused with a His tag at the C-terminus.

Bio-Activity

Shipping

Lyophilized proteins are shipped at ambient temperature.
Liquid proteins are shipped at frozen temperature with dry ice.
Upon receipt, store it immediately at the temperature recommended below.

Storage

Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.
After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.
Avoid repeated freeze/thaw cycles.

Operational Notes

For your safety and health, please wear a lab coat and disposable gloves for handling.

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn