



**Catalog No.: RP00976** **Recombinant**

Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.  
After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.  
Avoid repeated freeze/thaw cycles.


Validation Data



Recombinant Human IL-7RA/CD127 Protein was determined by SDS-PAGE under reducing conditions with Coomassie Blue.



Immobilized Recombinant Human IL7 at 2 µg/mL (100 µL/well) can bind Recombinant Human IL7R alpha, the EC<sub>50</sub> of Recombinant Human IL7R alpha is 48.44 ng/mL.



Immobilized CD127 Mouse mAb at 1µg/mL (25 µL/well) can bind Human CD34 with a linear range of 0.46-5.89ng/mL.