

Recombinant Mouse CD28 Protein

Catalog No.: RP01211 Recombinant

Sequence Information

Species Gene ID Swiss Prot Mouse 12487 P31041

Tags

C-hFc&His

Synonyms

CD28;Tp44;CD28

Product Information

Source

Purification

HEK293 cells

≥ 95 % as determined by SDS-

PAGE.

Calculated MW Observed MW

41.84 kDa 55-70 kDa

Endotoxin

< 0.1 EU/ μ g of the protein by LAL method.

Formulation

Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4.Contact us for customized product form or formulation.

Reconstitution

Centrifuge the vial before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stablizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

Contact

6		400-999-6126
\bowtie		cn.market@abclonal.com.cn
\odot	T	www.abclonal.com.cn

Background

Basic Information

Description

Recombinant Mouse CD28 Protein is produced by HEK293 cells expression system. The target protein is expressed with sequence (Asn20-Lys149) of mouse CD28 (Accession #NP_031668.3) fused with a Fc, 6×His tag at the C-terminus.

Bio-Activity

1.Measured by its binding ability in a functional ELISA. Immobilized Mouse CD86 Protein at 5 μ g/mL (100 μ L/well) can bind Mouse CD28 with a linear range of 0.0012-1.154 μ g/mL.

Shipping

The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below.

Storage

Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.

After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.


Avoid repeated freeze/thaw cycles.

Operational Notes


For your safety and health, please wear a lab coat and disposable gloves for handling.

^{*} ODODODODODODODO

Validation Data



Recombinant Mouse CD28 Protein was determined by SDS-PAGE under reducing (R) and non-reducing (NR) conditions.



Immobilized Mouse CD86 Protein at 5 μ g/mL (100 μ L/well) can bind Mouse CD28 with a linear range of 0.0012-1.154 μ g/mL.