




Validation Data



Recombinant Human Glutathione S-transferase omega-1/GSTO-1/SPG-R Protein was determined by SDS-PAGE under reducing conditions with Coomassie Blue.



Immobilized recombinant Human GSTO1 at 1 μ g/mL (100 μ L/well) can bind GSTO1 Rabbit mAb with a linear range of 0.5-4.33ng/mL.



Immobilized recombinant Human Annexin A5 at 10 μ g/mL (100 μ L/well) can bind Human GSTO1 with a linear range of 0.15-3.16 μ g/mL.