

Catalog No.: RP01569 **Recombinant**

Species	Gene ID	Swiss Prot
Drosophila melanogaste	31808	Q9W3C4

C-His

CG1787; Dmel\CG1787; DmHex2; HEX 2;
HEX2; HEXO2;Hexo2

Source	Purification
HEK293 cells	≥ 95 % as determined by SDS-PAGE.

< 0.1 EU/μg of the protein by LAL method.

Lyophilized from a 0.22 μm filtered solution of PBS, pH 7.4.

Centrifuge the vial before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn


Recombinant *Drosophila melanogaster* Beta-N-acetylhexosaminidase/Hexo2 Protein is produced by HEK293 cells expression system. The target protein is expressed with sequence (Asp37-Leu622) of *drosophila melanogaster* Beta-N-acetylhexosaminidase/Hexo2 (Accession #NP_525081.1) fused with a 6×His tag at the C-terminus.

Lyophilized proteins are shipped at ambient temperature.
Liquid proteins are shipped at frozen temperature with dry ice.
Upon receipt, store it immediately at the temperature recommended below.

Store at -20°C
Avoid repeated freeze/thaw cycles.

For your safety and health, please wear a lab coat and disposable gloves for handling.

Validation Data



Recombinant *Drosophila melanogaster* Beta-N-acetylhexosaminidase/Hexo2 Protein was determined by SDS-PAGE under reducing conditions with Coomassie Blue.