

**Catalog No.: RP01576LQ** **Recombinant**

| Species | Gene ID   | Swiss Prot |
|---------|-----------|------------|
| Mouse   | 100039150 | A9R9W0     |

## C-His

Major urinary protein 15; Mup13;  
Gm2068;MUP15;Major urinary protein  
15; Mup13; Gm2068;MUP15

|               |                                   |
|---------------|-----------------------------------|
| <b>Source</b> | <b>Purification</b>               |
| HEK293 cells  | ≥ 95 % as determined by SDS-PAGE. |

< 0.1 EU/μg of the protein by LAL method.

Supplied as a 0.22 µm filtered solution in PBS, pH 7.4.

Recombinant Mouse MUP-15 Protein is produced by HEK293 cells expression system. The target protein is expressed with sequence (Glu19-Glu180) of mouse MUP-15 (Accession #NP\_001186933.1) fused with a 6×His tag at the C-terminus.

Lyophilized proteins are shipped at ambient temperature.  
Liquid proteins are shipped at frozen temperature with dry ice.  
Upon receipt, store it immediately at the temperature recommended below.

Store at -80°C  
Avoid repeated freeze/thaw cycles.

For your safety and health, please wear a lab coat and disposable gloves for handling.


 | 400-999-6126

✉ | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Validation Data

---



Recombinant Mouse MUP-15 Protein was determined by SDS-PAGE under reducing conditions with Coomassie Blue.