

Catalog No.: RP01948LQ **Recombinant**

Species	Gene ID	Swiss Prot
African swine fever virus		P34204

N-His-TrxA-sumo

Product Information

E. coli

≥ 85 % as determined by SDS-PAGE.

[illegible]

Please contact us for more information.

Supplied as 0.22 μ m filtered solution in PBS, 300 mM NaCl, 10% glycerol, pH 7.4

Reconstitution

Basic Information

Recombinant African swine fever virus P30 Protein is produced by *E. coli* expression system. The target protein is expressed with sequence (Ile86-Lys194) of African swine fever virus P30 (Accession #P34204) fused with no tag.

Shipping

Lyophilized proteins are shipped at ambient temperature.
Liquid proteins are shipped at frozen temperature with dry ice.
Upon receipt, store it immediately at the temperature recommended below.

Storage

Store at -80°C
Avoid repeated freeze/thaw cycles.

Operational Notes

For your safety and health, please wear a lab coat and disposable gloves for handling.


Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Validation Data



Recombinant African swine fever virus P30
Protein was determined by SDS-PAGE under
reducing conditions with Coomassie Blue.