




Validation Data



Biotinylated Recombinant Human DLL3 Protein was determined by SDS-PAGE under reducing conditions with Coomassie Blue.



The purity of Biotinylated Human DLL3 is greater than 90% as determined by SEC-HPLC.



Immobilized Biotinylated Human DLL3, His Tag at 0.5µg/ml (100µl/well) on the streptavidin precoated plate (5µg/ml). Dose response curve for Anti-DLL3 Antibody, hFc Tag with the EC50 of 3.1ng/ml determined by ELISA (QC Test).