



Validation Data



Recombinant Mouse TNFSF11/RANKL/CD254 Protein was determined by SDS-PAGE under reducing conditions with Coomassie Blue.



Mouse bone marrow-derived BMDM P3 were used. Cells were induced with M-CSF (RP01216SC, 50 ng/mL) for 3 days, followed by induction with M-CSF (RP01216SC, 50 ng/mL) plus RANKL (RP02134C, 100 ng/mL) for 5 days (95 h). (top) RAW264.8 cells were treated with RANKL (RP02134C, 100 ng/mL) and induced for 4 days (76 h). (below).