

# Biotinylated Recombinant Human Angiopoietin-like 3/ANGPTL3 (17-219) Protein

Catalog No.: RP02273B **Recombinant**

## Sequence Information

Species	Gene ID	Swiss Prot
Human	27329	Q9Y5C1

### Tags

C-His&Avi

### Synonyms

ANG-5; ANGPT5; ANL3;  
FHBL2;ANGPTL3;ANGPT5;ANL3;FHBL2

## Product Information

Source	Purification
HEK293 cells	≥ 95 % as determined by Tris-Bis PAGE.

Calculated MW	Observed MW
26.6 kDa	30-38 kDa

### Endotoxin

< 1 EU/μg of the protein by LAL method.

### Formulation

Lyophilized from a 0.22 μm filtered solution of PBS, pH 7.4. Normally 8% trehalose is added as protectant before lyophilization. Contact us for customized product form or formulation.

### Reconstitution

Centrifuge the tube before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

## Contact

 | 400-999-6126

 | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

## Background

ANGPTL3 is a secreted glycoprotein that is structurally related to the angiopoietins. Mature human ANGPTL3 contains an N-terminal coiled coil domain and a C-terminal fibrinogen-like domain. ANGPTL3 is expressed in the liver from early in development through adulthood. Acts in part as a hepatokine that is involved in regulation of lipid and glucose metabolism. Proposed to play a role in the trafficking of energy substrates to either storage or oxidative tissues in response to food intake.

## Basic Information

### Description

Biotinylated Recombinant Human Angiopoietin-like 3/ANGPTL3(17-219) Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Ser17-Lys219) of Human Angiopoietin-like 3/ANGPTL3 fused with His tag and Avi tag at the C-terminus.

### Bio-Activity

### Shipping

The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below.

### Storage

Store at -20°C. Store the lyophilized protein at -20°C to -80°C up to 1 year from the date of receipt.

After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.


Avoid repeated freeze/thaw cycles.

### Operational Notes


For your safety and health, please wear a lab coat and disposable gloves for handling.



## Validation Data



Biotinylated Recombinant Human Angiotensin-like 3/ANGPTL3 (17-219) Protein was determined by Tris-Bis PAGE under reducing conditions.



Immobilized Anti-ANGPTL3 Antibody, hFc Tag at 1µg/ml (100µl/well) on the plate. Dose response curve for Biotinylated Human ANGPTL3, His Tag with the EC<sub>50</sub> of 5.0ng/ml determined by ELISA.