




Validation Data



Recombinant Cynomolgus/Rhesus macaque TNFRSF7/CD27 Protein was determined by Tris-Bis PAGE under reducing conditions.



The purity of Cynomolgus/Rhesus macaque CD27 is greater than 95% as determined by SEC-HPLC.



Immobilized Cynomolgus CD27 Ligand, His Tag at 2µg/ml (100µl/well) on the plate. Dose response curve for Cynomolgus/Rhesus macaque CD27, hFc Tag with the EC₅₀ of 41.5ng/ml determined by ELISA.